

Indikationer og protokoller

Indikation	Unit	FOV (cm)	Opløsning	Ikon	Dosis
Kæbeled	X1	5x15	0,200 (std)		57 μSv
Retinerede tænder og resorptioner i relation hertil	X1	4x4	0,150 (high)		42 μSv
Implantater – OK	X1	4x4 / 5x5	0,200 (std)		26 - 36 μSv
Implantater – UK	X1	4x4 / 5x5	0,150 (high)		42 - 60 μSv
Cyster/tumorer	X1	4x4 / 5x5	0,200 (std)		26 - 36 μSv
Endo	X1	4x4	0,075 (ultra)		56 μSv
Resorptioner (intern og cervikal)	X1	4x4	0,075 (ultra)		56 μSv
Orto-anomali	NewTom 5G	16x18	0,30 (reg)		190 μSv